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GH103 Flywheel Energy Storage Experimental System Product Specification

Edge Magnetic Pole Drive Technology | Flywheel Energy Storage Teaching & Research Platform

□ Product Overview

The GH103 flywheel energy storage experimental system adopts a W→A (Wheel-to-Arm) edge drive architecture, which differs from the traditional shaft-center drive mode. Electromagnetic force acts on the outer edge of the flywheel, and the moment arm is increased by 4-8 times compared to shaft-center drive (depending on flywheel diameter). It can be used to simulate the outer-rotor motor drive characteristics of large-scale energy storage flywheels. Suitable for university teaching demonstrations, scientific research experiments, and prototype validation..

□ Core Technology

● Drive Mode Comparison

Comparison Item	Traditional Shaft-Center Drive	GH103 Edge Symmetric Drive
Mechanical Model	$T = F \times r$	$T = F_t \times 2R$
Lever Arm	r (shaft radius, small)	2R (flywheel diameter, large)
Lever Arm Length	~10-20 mm	200 mm
Torque Enhancement	Baseline (1×)	4-8×

● Technical Features

- Proprietary Edge Drive: Electromagnetic force acts on flywheel outer rim, lever arm increased 4-8×
- Symmetric Couple Arrangement: 8 poles uniformly distributed, forming 4 symmetric couples
- Authentic Simulation: Perfectly reproduces outer-rotor motor characteristics of large-scale energy storage flywheels

□ Theoretical Formula System

● Basic Dynamic Equations

Formula Name	Expression	Physical Meaning
Law of Rotation	$\tau = J\alpha$	Torque equals moment of inertia times angular acceleration
Torque Definition	$T = F \times R$	Torque equals force times lever arm
GH103 Edge Torque	$T = F_t \times 2R$	Tangential force times diameter (couple arm), torque enhanced 4-8×

• Electromagnetic Drive Equations

Formula Name	Expression	Physical Meaning
Total Tangential Electromagnetic Force	$F_t = \eta \times f \times n \times \alpha$	Combined air gap, suction force, pole count and work coefficient for tangential drive force
Air Gap Permeance Coefficient	$\eta = 1/k$	Reciprocal of air gap distance, determines magnetic coupling strength
Output Torque	$T = F_t \times 2R$	GH103 core drive formula

• Kinematic Equations

Formula Name	Expression	Physical Meaning
Speed-Frequency Relation	$n = 60f / N_{poles} = 60f / 8$	Converts pulse frequency to flywheel speed (RPM)
Angular Velocity	$\omega = 2\pi n / 60 = \pi f / 4$	Converts speed to angular velocity (rad/s)
Stored Kinetic Energy	$E_k = \frac{1}{2} J \omega^2 = \frac{1}{2} J (2\pi n / 60)^2$	Kinetic energy stored in flywheel rotation (J)

• Energy Conversion Equations

Formula Name	Expression	Physical Meaning
Input Electrical Energy	$E_{elec} = \int_0^T V(t)I(t)dt$	Total electrical energy consumed by drive system
Conversion Efficiency	$\eta = E_{kinetic} / E_{elec} \times 100\%$	Electrical-to-kinetic energy conversion efficiency
Power Relation	$P = dE_k/dt = \tau \cdot \omega$	Instantaneous power equals torque times angular velocity

□ System Parameters

• Electromagnetic Drive System

Parameter	Value	Description
Operating Voltage	DC 12-100V	Rated operating voltage
Single Electromagnet Suction Force	0.5 kgf (4.9 N)	@DC12V, 25mm×20mm
Electromagnets per Pole	2 pcs	Symmetric arrangement
Total Pole Count	8 poles	Uniformly distributed, 45° interval
Electromagnetic Air Gap	0.8 mm	Gap between electromagnet and flywheel outer rim
Air Gap Permeance Coefficient	1.25 mm^{-1}	$\eta = 1/k$

● Mechanical Structure Parameters

Parameter	Value	Description
Flywheel Radius R	100 mm	Lever arm reference
Couple Arm 2R	200 mm	Diameter, torque action arm
Electromagnet Dimensions	25mm×20mm	Single electromagnet outline dimensions

□ Electromagnetic Force Calculation

● Total Tangential Electromagnetic Force

$$F_t = \eta \times f \times n \times \alpha$$

● Parameter Definitions

Symbol	Value	Physical Meaning	Unit
η	1.25	Air gap permeance coefficient, $\eta = 1/k$	mm^{-1}
k	0.8	Electromagnetic air gap distance	mm
f	9.8	Electromagnetic suction force per pole (equivalent to 1 kgf, 2×0.5 kgf)	N
n	8	Total number of magnetic poles	pcs
nm	2	Number of electromagnets per pole	pcs
α	0.5	Effective work coefficient ($p/2$, only half poles work tangentially)	—

● (@DC12V) / Calculation Example (@DC12V)

Total Tangential Electromagnetic Force Calculation:

$$\begin{aligned} F_t &= \eta \times f \times n \times \alpha \\ &= 1.25 \times 9.8 \times 8 \times 0.5 \\ &= 1.25 \times 9.8 \times 4 \\ &= 49 \text{ N} \end{aligned}$$

Output Torque Calculation: :

$$\begin{aligned} T &= F_t \times 2R \\ &= 49 \text{ N} \times 0.2 \text{ m} \\ &= 9.8 \text{ N}\cdot\text{m} \end{aligned}$$

Note: 8 poles uniformly distributed (45° per pole), 2 electromagnets per pole in symmetric arrangement, 16 electromagnets total forming 4 symmetric couples.

□ Experimental Formula Quick Reference Card

GH103 Experimental Formula Quick Reference Card		
Torque Calculation	$T = F_t \times 2R = 9.8 \text{ N}\cdot\text{m}$	
Inertia Calculation	$J = \frac{1}{2}mR^2$ (Disk Model)	
Kinetic Energy Calculation	$E_k = \frac{1}{2}J\omega^2$ [J]	
Speed Conversion	$\omega = 2\pi n/60 = 0.1047n$ [rad/s]	
Efficiency Calculation	$\eta = E_k / \int VI dt \times 100\%$	
PWM Duty Cycle	$D = t^+ / (t^+ + t^-)$	
Step Speed	$n = 60f/8$ [RPM]	

□ Eight-Pole Sequential Pulse Drive

- 8 magnetic poles uniformly distributed circumferentially, stepwise accelerating rotating magnetic field
- 8 steps per revolution, step angle 45°, smooth operation with low vibration
- Frequency 10-400Hz stepless adjustment, precise speed control

□ Technical Parameters

● Mechanical System

Parameter	Specification	Description
Flywheel Diameter	φ200 mm	High-strength aluminum alloy, dynamic balance G2.5
Weight Code System	5kg + 3kg + 2kg	Modular combination, inertia 0.025-0.25 kg·m ² adjustable
Pole Configuration	Single stage 8 poles	Voltage-adjustable electromagnetic coil, φ25mm
Lever Action Angle	45° (π/4)	Optimized torque transfer efficiency
Air Gap Adjustment	0.2-2 mm	Adaptable to different experimental requirements
Overall Dimensions	600 × 600 × 400 mm	Laboratory desktop design
Overall Weight	~380 kg	Stable anti-vibration base

● Electrical Drive

Parameter	Specification	Description
Drive Voltage	DC 12-100 V	Wide range adjustable, adaptable to multiple power sources
Drive Voltage	10-400 Hz	Corresponding speed 75-3000 RPM ($n = 60f/8$)
Positive Pulse Width	1-20 ms	Energy injection time adjustable
Negative Pulse Width	1 ms	Fixed freewheeling/braking
Peak Power	1.5 kW	Meets acceleration and steady-state requirements
Input Power	AC 100-240V	Wide voltage, global universal

● Motion Performance

Parameter	Specification	Description
Speed Range	10-3000 RPM	Covers low-speed debugging to high-speed energy storage
Rated Speed	600 RPM	Standard teaching condition ($f = 80\text{Hz}$)
Maximum Linear Velocity	31.4 m/s	Edge linear velocity at 3000 RPM ($v = \omega R$)
Angular Velocity	1.05-314 rad/s	Wide range dynamic response
Angular Acceleration	0.02-2 rad/s ²	Precision acceleration control ($\alpha = \tau/J$)
Start-up Time	<8s (0-600 RPM)	Rapid experimental condition establishment


● Measurement & Control System

Parameter	Specification	Description
Speed Detection	Hall sensor + magnetic cone	8 pulses/revolution, non-contact measurement
Waveform Analysis	Tektronix TDS7404B	100MHz, 4-channel real-time observation
Electrical Measurement	HAMEG HM1507-3	Voltage/current/frequency/duty cycle
Speed Display	4-digit LED display	Intuitive RPM reading
Safety Protection	GH19C Intelligent System	Five-fold protection: overspeed/overcurrent/overheat/vibration/displacement

□ **Ten-Dimensional Experimental Variables**

No.	Variable Name	Symbol	Adjustment Range	Experimental Significance	Related Formula
1	Counterweight Mass	m	2–10 kg	Moment of inertia J adjustment	$J = \sum m_i r_i^2$
2	Electromagnetic Force	F	0–500 N	Drive torque $T = F \cdot R$	$T = F \cdot R$
3	Drive Voltage	V	12–100 V	Control excitation current	$V = IR + L(di/dt)$
4	Drive Current	I	0–15 A	Electromagnetic force excitation source	$F \propto B \propto I$
5	Pulse Frequency	f	10–400 Hz	Speed setting	$n = 60f/8$
6	Positive Pulse Width	t ⁺	1–20 ms	Energy injection duration	$E_{pulse} = VI t^+$
7	Angular Acceleration	α	Measured value	Dynamic response output	$\alpha = \tau/J$
8	Magnetic Field Displacement Angle	θ	±7–13.5°	Optimal action angle	$F_t = F \cdot \sin\theta$
9	Magnetic Induction Intensity	B	0–1.2 T	Proportional to current	$B = \mu_0 \mu_r NI / l_0$
10	Flywheel Phase	(x,y)	$[0, 2\pi]^2$	State space variable	$\dot{x} = Ax + Bu$

□ **Core Experimental Functions**

● **Basic Physics Experiments**

- Law of Rotation Verification: $\tau = J\alpha$, linear torque-angular acceleration relationship
- Inertia Measurement Experiment: $J = \sum m_i r_i^2$, compare theoretical calculation with measured values
- Energy Conservation Demonstration: $E_{elec} \rightarrow E_k \rightarrow E_{elec}$ conversion efficiency analysis

● **Motor Drive Experiments**

- Switched Reluctance Motor Principle: 8-pole sequential in-phase mechanism, $n = 60f/8$
- Mechanical Characteristic Curves: Speed-torque-efficiency three-dimensional relationship
- PWM Control Technology: $D = t^+ / (t^+ + t^-)$, pulse width modulation and duty cycle optimization

● **Energy Storage System Experiments**

- Charge-Discharge Cycles: $\eta = E_k / \int VI dt$, energy storage and release efficiency analysis
- Dynamic Response Testing: Step/ramp/disturbance response, $\alpha = d\omega/dt$
- Braking Energy Recovery: Comparison of dissipative braking and plug braking

● Advanced Research Experiments

- Multivariable Coupling Analysis: Ten-dimensional parameter orthogonal experimental design
- Edge Drive Optimization: Air gap/angle/frequency collaborative optimization
- Model Identification and Validation: Theoretical mathematical modeling vs. experimental comparison

□ Target Users

User Type	Application Scenario	Value Proposition
Universities	Physics/Electrical/Energy major experimental teaching	Fills the gap in flywheel energy storage experimental equipment, integrating theory with practice
Vocational Colleges	New energy technology practical training	Master energy storage system installation, commissioning and maintenance skills
Research Institutes	Flywheel energy storage mechanism and control research	Provides a repeatable precision experimental platform
Enterprise R&D	Motor drive and energy management pre-research	Low-cost validation of innovative solution feasibility

□ Product Advantages

Comparison Item	raditional Experimental Equipment	GH103
Drive Mode	Shaft-center motor drive	Edge magnetic pole drive, $T = F_t \times 2R$, torque enhanced 4-8×
Parameter Adjustment	Fixed parameters, single experiment	Ten-dimensional adjustable, flexible experimental design
Measurement Means	Mechanical/simple electronic	Hall sensor + digital oscilloscope, precision acquisition
Safety Protection	Basic protection	GH19C five-fold intelligent protection, safe and reliable
Expansion Capability	Closed structure	Modular design, supports multi-level expansion and vacuum upgrade
Theoretical Support	No systematic formula system	Complete formula system, integrated teaching and research

□ Standard Configuration

- ✓ GH103 Main Unit (flywheel rotor, 8-pole drive, control system, protective cover)
- ✓ Counterweight Set (flywheel weights: 5kg×1, 3kg×1, 2kg×1, with quick-lock mechanism)
- ✓ Tektronix TDS7404B Digital Oscilloscope (100MHz, 4-channel)
- ✓ HAMEG HM1507-3 Multifunction Measuring Instrument
- ✓ Hall Sensor Kit (8 sets, including magnetic cones)
- ✓ Electrical Connection Harness (power lines, signal lines, BNC cables)
- ✓ Operation Manual (including experimental guidance and theoretical derivation)
- ✓ Factory Inspection Certificate and Warranty Card

□ Optional Extensions

Option	Function Description	Application Scenario
Vacuum Chamber Kit	≤100Pa / Stainless steel vacuum cover + vacuum pump, vacuum ≤100Pa	High-speed low wind resistance experiments, wind loss mechanism research
Secondary Expansion Module	Add 8 magnetic poles, 16-pole configuration	High-resolution step drive, micro-vibration analysis
Data Acquisition System	16-bit USB-DAQ, 4-channel synchronization, supporting software	Automated experimental recording, large data volume acquisition
Temperature Monitoring Module	8-channel PT100 temperature sensors, magnetic pole thermal management	Long-term operation thermal characteristics research, efficiency optimization
Ceramic Bearing Upgrade	High-speed bearings, supports 5000 RPM expansion	Ultra-high-speed energy storage experiments, limit performance testing

□ Service Commitment

Service Item	Content
Warranty Period	1 year for complete machine, 3 years for controller and driver
Installation & Commissioning	Free on-site installation, commissioning, and operation training
Technical Support	7×24 hour telephone support, remote diagnosis, on-site service
Software Upgrade	Continuous optimization of control algorithms, lifetime free upgrades
Custom Development	Customize parameters according to special requirements, develop exclusive experimental solutions

□ Technical Specification Quick Reference

Key Indicator	Value
Flywheel Diameter	200 mm
Speed Range	10-3000 RPM
Drive Voltage	12-100 V
Peak Power	1.5 kW
Moment of Inertia	0.025-0.25 kg·m ²
Angular Acceleration	0.02-2 rad/s ²
Overall Weight	~380 kg
Protection Grade	IP20

□ Contact Us

Company: GeHa (Shenzhen) Technology Co., Ltd

Address: No. 1301, Guanlan Street, Longhua District, Shenzhen, China

Telephone: +86-755-28168358

Fax: +86-755-28168293

Website: <http://www.szgeha.com>

Email: szgeha@163.com

GH103 Flywheel Energy Storage Experimental System — Transforming abstract energy storage theory into intuitive experimental data. Edge magnetic pole drive, empowering energy storage teaching practice.